

www.imv-technologies.com

Next Generation CASA System



R



Next Generation CASA Systems IVOS II and CEROS II

- Standardization
- Objective & reliable measurement
- Accuracy and repeatability
- Traceability
- Multiparameters
- Decision-making tool

Software Overview

Interactive Illumination Check



The new software features the unique color-coded Illumination Check to optimize identification of the sperm head and sperm tail. This feature removes any guesswork in setting the illumination and promotes consistency between all users.

Screen Layout



- 1. Controls to initiate analysis
- 2. Quick selection of analysis setup
- 3. Results, updated in real time.
- 4. Live image area
- 5. Thumbnail gallery of captured video images.

Dose Input		₽.		
Sperm / Dose	3000.0			•
Dose Volume	80		(=)(•
Adjustments / I	Results			
Concentration	396.3			M
Total Count	85			
C Motile	92.9			
Progressive	69.4			1
Morph Normal			95.3	6
Processing			2 Aut	•
Extender Volume		2722		ml
Final Volume		3022		mi
Adjusted Concentration		39.4		Min
Number Of Doses		37		



Key benefits

- Automated morphology
- Dose calculation:
 Concentration
 - Extender volume Number of doses
- Automated heated stage *
- Fluorescence *
- Redesigned graphical user interface
- Intuitive navigation
- Multilanguage:

English, Russian, Chinese and Spanish (with additional languages to come)

High-speed digital camera:

Superior image quality Seamless image capture and playback

- Unlimited storage of predefined analysis setups
- Built-in database with customizable reports
- ASCII data export
- Fast: 8 seconds / 10 fields (30 images per field)

*Available with IVOS II



Analysis Results

Real time Updating

Analysis results and dose calculation are updated in real time as each field is added to the analysis.

Full Field Playback

User-definable. Selecting a thumbnail image from the gallery opens the captured video. The video may be replayed in full or you may scroll through, frame by frame.

Zoom Cell Playback

Selecting a track from the playback image opens up a zoomed image of the tracked cell. Kinematic measurements are displayed relative only to the selected track.

Cell Morphometry and Automated Morphology

The software also provides cell morphometry values, including head length, width, perimeter and area. It also identifies automatically:

- Distal droplets
- Proximal droplets
- Distal Midpiece Reflexes
- Coiled tails
- Bent tails



Measured and Calculated Parameters

Counts, concentration and percentages

Total, static, motile, progressive and slow cells

Doses

- Extender volume
- Final volume
- Number of doses

Kinematic measurements

DAP, DSL, DCL, VAP, VSL, VCL, ALH, STR, LIN, BCF, WOB (Average, SD, Median)

Morphometric characteristics

Head length, width, perimeter and area. Tail length and straightness. Droplet distance.

Automated morphology

- Distal droplets
- Proximal droplets
- Distal Midpiece Reflexes
- Coiled tails
- Bent tails





Disposable Leja counting chambers offer you the most reliable, quick, and accurate solution!

- Best brand name on the market
- Specially designed for semen analysis
- Tested not spermotoxic
- Disposable slides: time-saving
- Highest precision, repeatability and reliability
- Accurate and constant chamber depth (20 micron)
- Easy-to-use
- Available with 4 or 8 chambers



Three quality control stages per batch:

- Chamber height* Visual inspection (18 points) Toxicity control
- * Certificate of analysis for each batch available at www.lejaslides.com





Semen analysis by IMV Technologies



IVOS II

Ref. 024911: IVOS II Ref. 024910: IVOS II fluorescence

- Built-in system
- Automated heated stage
- Fluorescence
- Automated morphology

ISO 9001:2008 and medical reference ISO 13455:2003 certified Production site complies with the strictest standards of quality and safety Continuous improvement of product, production and quality control processes All materials and finished products potentially in contact with living cells tested for bio compatibility 95% made in France



CEROS II

Ref. 024905: CEROS II Zeiss Ref. 024904: CEROS II Olympus

- Manual heated stage
- Automated morphology







EasyCyte II plus Ref.022841



EasyCyte II mini Ref.022843



www.imv-technologies.com

+ 33 (0) 233 346 464 contact@imv-technologies.com