



www.imv-technologies.com

GT DEEP™ concept

Post-cervical Insemination by IMV




GT DEEP™ Technology

- Proven since 1999
- Maximises genetic potential
- Easy-to-use
- A complete concept, from automated packaging to insemination
- Full traceability




GT DEEP™ concept


Post-cervical Insemination by IMV



Beginning of Artificial Insemination (AI)



1 minute after AI





2 minutes after AI

GT DEEP technology

A pioneer of post-cervical artificial insemination (PCAI) of sows, IMV Technologies has developed a full concept of post-cervical insemination, from packaging to insemination. This technique, patented in 1985, allows the introduction of semen 20 cm beyond the cervix, and limits spermatozoa loss.

Thanks to this experience, know-how and many partnerships, IMV Technologies' PCAI concept has become a worldwide standard.



GTB Bags

Perfect fit with any type of catheter

- Direct connection without additional connector
- Leak-proof connection

Perfect preservation conditions

- Two-layered film with guaranteed non-spermicidal properties
- Increases conservation time (flat storage)
- Sanitary: tip is protected until AI

Compact and practical

- Space saving of up to 15%
- Easy-opening, no tool required

GTB Bags ¹

Capacity: up to 90 ml
Design optimized for traditional insemination and hands-free use

Mini GTB Bags ²


Capacity: up to 40 ml
Special design for PCAI
Allows dose emptying by simple pressure



IMV Catheters

As a leader in this field, IMV Technologies has developed a range of PCAI catheters that combines ease-of use, efficiency and quality of manufacturing.

- Atraumatic: PCAI catheter tip is completely smooth, without any asperities
- Sperm friendly: tip assembly without glue
- White tip to display the presence of any blood
- Available with or without connector



Automated Packaging

IMV offers a full line of packaging systems for GTB bags.

- Flexible: fully compatible with GTB or Mini GTB bags

- Guaranteed traceability
- Sanitary: easy to clean, and optimal use of disposable solutions
- Easy-to-use and precise (+/-1ml)
- Reliable

GTB 1000™ ¹

Automated filling, sealing and labeling machine

GTB 600™ ²

Semi-automated filling and sealing

GTB 250™ ³

Manual filling and sealing



Bibliography

P.F. Watson, J.R. Behan – Intrauterine insemination of sows with reduced sperm numbers: results of a commercially based field trial – *Theriogenology* 57 (2002) 1683-1693

P.F. Watson, J.R. Behan – A field investigation of intra-cervical insemination with reduced sperm numbers in gilts – *Theriogenology* 66 (2006) 338-343

Sasha Gibson, MSc; Robert J. Tempelman, PhD; Roy N. Kirkwood, DVM, PhD – Effect of oxytocin supplemented semen on fertility of sows bred by intrauterine insemination – *Journal of Swine Health and Production* – July and August, 2004

Donald G. Levis, Scott Burroughs and Sara Williams – Use of Intra-Uterine Insemination of Pigs: Pros, Cons & Economics – *The Ohio State University, Columbus, Ohio 43210-1095, Danbred USA, Dorchester, Nebraska 68343, Institute of Theriogenology, Faculty of Veterinary Science, National University of La Plata, La Plata, Argentina*

Kakanang Buranaamnuy, Termpong Wongtawan, Sutthapit Masuwatana, Padet Tummaruk, Mongkol Techakumphu
Intra-uterine and Deep intra-uterine Insemination using Cryopreserved Boar Semen in Spontaneously-ovulating Sows
Thai J. Vet. Med. 2010 40(2): 215-219.

Results from P.F. Watson & J.R. Behan 2002

Table 1

A comparison of fertility and fecundity of sows after artificial insemination twice in a single oestrus with a standard device (Goldenpig®: control) and a novel deep insemination device (Deep Goldenpig™)

Device	Sperm dose (Billions)	Pregnancy (%)	Farrowing (%)	Litter size (n)	Number born live (n)
Goldenpig®	1	66.2	65.8	10.3	9.0
	2	91.1	91.8	12.6	10.9
	3	91.3	91.1	12.5	10.9
Deep Goldenpig™	1	88.7	86.9	12.1	10.9
	2	92.6	92.5	12.3	10.8
	3	91.8	90.5	12.3	11.0
Total (N)		3 230	3 201	2 768	2 768

Results from P.F. Watson & J.R. Behan 2006

Table 2

A comparison of fertility and fecundity of gilts after insemination twice in a single oestrus with the novel gilt catheter with different sperm cell numbers

Sperm dose (x10 ⁹)	Pregnant gilts (n)	Pregnancy (%)	Farrowing (%)	Mean Litter size (n)
1	553	92.3	90.9	11.39
2	559	93.2	91.1	11.28

Table 3

Reproductive performance of sows inseminated with Goldenpig® (3.0 billion sperm cells) or DeepGoldenpig™ (1.5 billion sperm cells)

Sperm per dose	Sows cycling by 7 days after weaning			Avg piglets born per litter		Fecundity index (born alive)
	Number of sows		Farrowing rate, %	Total	Born alive	
Goldenpig® (3.0 billion)	Inseminated	Farrowed				
	192	149	77.60	11.86	10.28	798
Deep Goldenpig™ (1.5 billion)	189	144	76.19	11.70	10.55	779
Difference	3	5	1.41	0.16	0.27	19




- ISO 9001:2008 and medical reference ISO 13485:2003 certified
- Production site complies with the strictest standards of quality and safety
- Continuous improvement of product, production and quality control processes
- All materials and finished products potentially in contact with living cells tested for bio compatibility
- 95% made in France



+ 33 (0) 233 346 464
contact@imv-technologies.com

