

Swine A.I. catheter Gold standard

- **Improved performance:**
 - Perfect lock into the cervix
 - Tight seal assures reduced backflow of semen
- **Reliable-foam tip permanently fixed:**
 - Foam tip is molded on to assure the highest success for the breeder
- **Proven technology:**
 - First molded "Foam Tip" catheter introduced in 1989
 - More than 400 million inseminations performed with IMV's Goldenpig
- **Design:**
 - Foam tip prevents trauma to the genital track
- **Bio-secure & consistent quality:**
 - Manufactured under strict sanitary protocol
 - Manufactured under consistent ISO 9001 procedures
- **Complete line:**
 - IMV offers the Goldenpig® for sow;
 - Goldengilt® for gilts;
 - and the DeepGoldenpig Concept III® for post cervical inseminations



Goldenpig®



Goldengilt®

Made in France
ISO 9001



DeepGoldenpig
Concept III®

Often imitated, never duplicated

Ask for the original Goldenpig®

WORLD LEADER IN REPRODUCTION BIOTECHNOLOGIES



www.imv-technologies.com

SWINE A.I. CATHETER GOLD STANDARD:


Made in France - ISO 9001

Reproduction is a key success factor for swine producers. **Your profitability depends highly on your herd's reproductive performance.**

IMV Technologies has been involved with swine AI for over 25 years, and we are aware of the importance of AI, and AI's impact on producer's bottom line. Our daily objective is to assure the supply of the highest quality products.

We can control our quality standards because **we manufacture all of our products in our plants** under strict quality controls, certified ISO 9001 procedures, with trained workers.

This is the reason why **IMV Technologies is today's best choice for supplier of AI catheters.**



The specific design of the Goldenpig® catheter has been developed to insure that the sow will not be harmed. The foam around the tip of the catheter protects the sow's cervix, and provides a perfect lock at the same time.

Goldenpig®

Since its introduction **more than 400 million inseminations** have been performed using the Goldenpig® catheter! The proven design of the product makes the Goldenpig® the ideal catheter for the reproductive tract of the sow or gilt. It is easy to insert and easy to use! No need to worry about it "locking-in" in the cervix. The foam tip works to decrease semen backflow without inhibiting uterine contractions. The size of the tip makes it too large to accidentally penetrate the urethra.

Goldengilt®

IMV offers a catheter designed specifically for gilts. The catheter has a smaller diameter, greater taper, and slightly longer foam tip. The Goldengilt® catheter is **ideal for the smaller reproductive tract of gilts!**

DeepGoldenpig Concept III®

The DeepGoldenpig® was created for successful "deep" or trans-cervical intra-uterine AI. The concentration of sperm cells per dose can be significantly reduced with this technique as the semen is placed deeper into the reproductive tract. Swine producers using the DeepGoldenpig® catheter are successful in using **a reduced sperm cell concentration that result in excellent conception rates.** A small flexible catheter is inserted through the conventional Goldenpig® which enables the user to navigate through the female's cervix. The small flexible inner catheter is tipped with a protective rounded head which removes the risk of uterine tract damage. Semen can then be deposited past the cervix into the body of the uterus or within the uterine horns.

WORLD LEADER IN REPRODUCTION BIOTECHNOLOGIES



www.imv-technologies.com