
## Improve your Conception Rate - IFT SmartWarmer™

Embryo/Semen/IVF/Sexed Semen/Oocytes






## Smart Over Heating Protection / Cut-off system.

No matter which temperature is selected, the microcontroller will automatically cease heating when the sensors detect a sudden increase of heat. LED lights will light up, alerting the user.





Nominal capacity:

Charging time:

Operating time:

10,400mAh (37Wh)

5°C ambient - set @ 36°C / 7hrs continuous

6hrs @ 2A









IFT Instruments SrI - info@iftinstruments.com - www.iftinstruments.com